

miProfile™ Human miRNome miRNA qPCR Array(384-well)

For focused group profiling of human miRNA whole genome expression

Cat. No. QM042-F (5 x 96-well plate, Format F)

Cat. No. QM042-G (5 x 96-well plate, Format G)

Cat. No. QM042-H (5 x 96-well plate, Format H)

Available as 1 set or 6 sets. Each set contains 1,700 unique miRNA primers deposited in five 384-well plates.

Introduction

The miProfile human miRNome miRNA PCR array (384-well) is a set of five 384-well plates, covering 1,700 of best characterized and annotated miRNAs based on miRBase V18. Each 384-well plate contains up to 360 pairs of PCR primers (forward: miRNA-specific primer; reverse: universal primer), which are pre-deposited in each well. Each plate also designated 24 control wells for monitoring the efficiency of every step of the experiment - from reverse transcription to qPCR reaction.

- QM042 plate 01: 360 unique miRNA PCR primer pairs
- QM042 plate 02: 360 unique miRNA PCR primer pairs
- QM042 plate 03: 360 unique miRNA PCR primer pairs
- QM042 plate 04: 360 unique miRNA PCR primer pairs
- QM042 plate 05: 260 unique miRNA PCR primer pairs

Shipping and storage condition

Shipped at room temperate

Stable for at least 6 months when stored at -20°C

Array format

GeneCopeia provides three qPCR array formats (F, G, and H) suitable for use with the following real-time cyclers.

Important note: Upon receiving, please check to make sure that the correct array format was ordered to ensure the compatibility with your qPCR instrument.

Plate format	Instrument provider	qPCR instrument model
F (384-well)	Applied Biosystems	7900HT (384-well block), ViiA™7 (384-well block)
G (384-well)	Bio-Rad Laboratories	CFX384™ (384-well block)
H (384-well)	Roche Applied Science	LightCycler® 480 (384-well block)

Quality control

1. Each miRNA-specific primer in the miProfile miRNA qPCR array has been experimentally validated to yield a single dissociation curve peak and to generate a single amplicon of the correct size for the targeted miRNA.
2. The positive PCR controls (PCR) have been verified to amplify a single amplicon of the correct size with Ct values around **20±2**.

3. The Spike-in reverse transcription controls (RT) have been verified to amplify a single amplicon of the correct size with Ct values around **20±3**.
4. $R_2 > 0.99$ was observed for high inter/ intra-array reproducibility.

Materials required but not provided

All-in-One™ miRNA First-Strand cDNA Synthesis Kit
 All-in-One™ qPCR Mix
 Total RNA extraction kit (RNAzol® RT RNA extraction reagent is recommended)
 DNase/RNase free tips, PCR reaction tubes, 1.5 ml microcentrifuge tubes
 5 ml and 10 ml graduated pipettes, beakers, flasks, and cylinders
 10 µl to 1,000 µl adjustable single channel micropipettes with disposable tips
 5 µl to 20 µl adjustable multichannel micropipette, disposable tips, and reservoir
 qPCR instrument, compatible with gene qPCR arrays ordered

Array layout

	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24
A	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24
B	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48
C	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64	65	66	67	68	69	70	71	72
D	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	88	89	90	91	92	93	94	95	96
E	97	98	99	100	101	102	103	104	105	106	107	108	109	110	111	112	113	114	115	116	117	118	119	120
F	121	122	123	124	125	126	127	128	129	130	131	132	133	134	135	136	137	138	139	140	141	142	143	144
G	145	146	147	148	149	150	151	152	153	154	155	156	157	158	159	160	161	162	163	164	165	166	167	168
H	169	170	171	172	173	174	175	176	177	178	179	180	181	182	183	184	185	186	187	188	189	190	191	192
I	193	194	195	196	197	198	199	200	201	202	203	204	205	206	207	208	209	210	211	212	213	214	215	216
J	217	218	219	220	221	222	223	224	225	226	227	228	229	230	231	232	233	234	235	236	237	238	239	240
K	241	242	243	244	245	246	247	248	249	250	251	252	253	254	255	256	257	258	259	260	261	262	263	264
L	265	266	267	268	269	270	271	272	273	274	275	276	277	278	279	280	281	282	283	284	285	286	287	288
M	289	290	291	292	293	294	295	296	297	298	299	300	301	302	303	304	305	306	307	308	309	310	311	312
N	313	314	315	316	317	318	319	320	321	322	323	324	325	326	327	328	329	330	331	332	333	334	335	336
O	337	338	339	340	341	342	343	344	345	346	347	348	349	350	351	352	353	354	355	356	357	358	359	360
P	NC	NC	NC	NC	HK1	HK1	HK2	HK2	HK3	HK3	HK4	HK4	HK5	HK5	HK6	HK6	RT	RT	RT	RT	PCR	PCR	PCR	PCR

Figure1. Illustration of miProfile miRNA qPCR array (384-well plate)

- **miRNA primer pairs:** Wells 1-360 are designated wells for pre-deposited miRNA primer pairs.
- **NC:** Negative controls, which only have the pre-deposited reverse universal primers.
- **HK1-6:** Twelve pre-deposited housekeeping snRNAs primer pairs, which can be used as endogenous positive controls as well as for array normalization.
- **RT:** Four replicates of spike-in reverse transcription controls, which can be used to monitor the efficiency of the RT reaction. These pre-deposited primer pairs specifically amplify the cDNA template reversed transcribed from the spike-in exogenous RNA in the sample.
- **PCR:** Four replicates of positive PCR controls, which are used to verify the PCR efficiency by amplifying the pre-deposited DNA template with its specific pre-deposited primer pairs.

miRNA primer list

The primer list can be downloaded from <http://www.genecopoeia.com/product/qpcr-arrays/mirna/mirnomia.php>.

Limited Use License

Following terms and conditions apply to use of miProfile™ Human miRNome miRNA qPCR Array (384-well) (the Product). If the terms and conditions are not acceptable, the Product in its entirety must be returned to GeneCopoeia within 5 calendar days. A limited End-User license is granted to the purchaser of the Product. The Product shall be used by the purchaser for internal research purposes only. The Product is expressly not designed, intended, or warranted for use in humans or for therapeutic or diagnostic use. The Product must not be resold, repackaged or modified for resale, or used to manufacture commercial products or deliver information obtained in service without prior written consent from GeneCopoeia. This Product should be used in accordance with the NIH guidelines developed for recombinant DNA and genetic research. Use of any part of the Product constitutes acceptance of the above terms.

Limited Warranty

GeneCopoeia warrants that the Product meets the specifications described in the accompanying Product Datasheet. If it is proven to the satisfaction of GeneCopoeia that the Product fails to meet these specifications, GeneCopoeia will replace the Product. In the event a replacement cannot be provided, GeneCopoeia will provide the purchaser with a refund. This limited warranty shall not extend to anyone other than the original purchaser of the Product. Notice of nonconforming products must be made to GeneCopoeia within 30 days of receipt of the Product. GeneCopoeia's liability is expressly limited to replacement of Product or a refund limited to the actual purchase price.

GeneCopoeia's liability does not extend to any damages arising from use or improper use of the Product, or losses associated with the use of additional materials or reagents. This limited warranty is the sole and exclusive warranty. GeneCopoeia does not provide any other warranties of any kind, expressed or implied, including the merchantability or fitness of the Product for a particular purpose.

GeneCopoeia is committed to providing our customers with high-quality products. If you should have any questions or concerns about any GeneCopoeia products, please contact us at 301-762-0888.

© 2016 GeneCopoeia, Inc.

GeneCopoeia, Inc.
9620 Medical Center Drive, Suite 101
Rockville, MD 20850
+1 (301) 762-0888
+1 (866) 360-9531
inquiry@genecopoeia.com