

Genome-CRISP™ CRISPR 人类 Innate Kinases & Ubiquitin Ligases sgRNA 混合文库（菌液）

由靶向人类内生激酶及泛素连接酶基因的 sgRNA 组成。

简介

Genome-CRISP™ CRISPR 人类 innate kinases & ubiquitin ligases sgRNA 文库靶向 239 个编码人类内生激酶及泛素连接酶的基因。这些 sgRNA 的设计都经过优化，以保证它们靶向目的基因所有剪接体序列的公共区域。sgRNA 文库能通过基因敲除实现基因组水平的功能缺失型筛选，是哺乳动物细胞系统基因分析、促进基因研究、基因组级别功能调查（如信号转导通路）以及药物研究（如靶标甄别和药物机制研究）的强大工具。

- Cat. No. L01-LS03-B1-1A（菌液，混合文库，分文库 1A）
- Cat. No. L01-LS03-B1-2A（菌液，混合文库，分文库 2A）
- Cat. No. L01-LS03-B1-1B（菌液，混合文库，分文库 1B）
- Cat. No. L01-LS03-B1-2B（菌液，混合文库，分文库 2B）


文库由 4 个分文库组成。分文库 1A 包含 118 个 sgRNA；分文库 2A 包含 118 个 sgRNA；分文库 1B 包含 118 个 sgRNA；分文库 2B 包含 121 个 sgRNA。每分文库为 1 管体积极约 1 mL 的甘油菌种。

运输及储存条件

干冰运输。-80°C 条件下能稳定储存至少 2 年。

载体信息

载体名称: pCRISPR-LvSG03
载体类型: 慢病毒
筛选标记: 嘌呤霉素
启动子: U6
载体大小: 8415 bp



克隆信息

测序引物: TTCTTGGGTAGTTTGACG

有限使用许可

Following terms and conditions apply to use of Genome-CRISP™ CRISPR Human Innate Kinases & Ubiquitin Ligases sgRNA Library (the Product). If the terms and conditions are not acceptable, the Product in its entirety must be returned to GeneCopoeia within 5 calendar days. A limited End-User license is granted to the purchaser of the Product. The Product shall be used by the purchaser for internal research purposes only. The Product is expressly not designed, intended, or warranted for use in humans or for therapeutic or diagnostic use. The Product must not be resold, repackaged or modified for resale, or used to manufacture commercial products or deliver information obtained in service without prior written consent from GeneCopoeia. This Product should be used in accordance with the NIH guidelines developed for recombinant DNA and genetic research. Use of any part of the Product constitutes acceptance of the above terms.

有限质保生命

GeneCopoeia warrants that the Product meets the specifications described in the accompanying Product Datasheet. If it is proven to the satisfaction of GeneCopoeia that the Product fails to meet these specifications, GeneCopoeia will replace the Product. In the event a replacement cannot be provided, GeneCopoeia will provide the purchaser with a refund. This limited warranty shall not extend to anyone other than the original purchaser of the Product. Notice of nonconforming products must be made to GeneCopoeia within 30 days of receipt of the Product. GeneCopoeia's liability is expressly limited to replacement of Product or a refund limited to the actual purchase price.

GeneCopoeia's liability does not extend to any damages arising from use or improper use of the Product, or losses associated with the use of additional materials or reagents. This limited warranty is the sole and exclusive warranty. GeneCopoeia does not provide any other warranties of any kind, expressed or implied, including the merchantability or fitness of the Product for a particular purpose.

GeneCopoeia is committed to providing our customers with high-quality products. If you should have any questions or concerns about any GeneCopoeia products, please contact us at 301-762-0888.

© 2016 GeneCopoeia, Inc.

广州高新技术产业开发区
广州科学城掬泉路3号
广州国际企业孵化器F区8楼
510663

电话: 4006-020-200; 020-32068595

传真: 020-32052877

网址: www.igenebio.com

邮箱: sales@igenebio